LIFT mee

-DERMA PEEL-

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User manual

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**LIFT
mee**
-DERMA PEEL-

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**LIFT
mee**
-DERMA PEEL-

Congratulation to your purchase of the New LIFTmee DERMA PEEL

LIFTmee Microderm is a vacuum assisted peeling system, which massages and peels your skin at the same time. The massage stimulates the blood circulation of the skin giving a firmer, smoother and revitalising appearance. Through the exfoliation, the dead skin cells sitting on the skin's outmost surface and blemishes are effectively removed. The soft, smooth vacuum stimulates the renewal process revealing the fresher younger skin cells below. Microdermabrasion also makes it easier for your skin products to penetrate into the deeper layers of the skin to give a youthful appearance. Already after the first use your skin can feel smoother and fresher looking.

IMPORTANT NOTICE

Danger

Never use the device or adaptor in a wet environment. Keep adaptor dry.

Warning

Only use the appliance with adaptor supplied and original parts. Recommended only for persons above 12 years. This device contains a magnet. People with medical equipment (cardiac pacemakers etc) should refrain from use.

Caution

The equipment may only be used for the described cosmetic actions.
Not to be used in combination with a mechanical or chemical peeling product.
Maximum use of device, twice a week.
To be used only on healthy skin.
Not to be used if you have skin infections, sun burn warts, moles, scars cuts, abrasions, tattoos, skin inflammation, eczema, extremely sensitive skin or permanent makeup.
Consult your doctor or dermatologist for advice.

General advice

For a minimum of three days after use apply sun protection to treated areas.
Do not use prior to intensive sun bathing.
Charge the device fully at least every 3 months to increase the life cycle of the batteries.

Preparation

The device should be charged for a minimum of 8hrs before use. The skin should be clean and dry. Before first use on your face, test skin reaction. Best on your forearm. Before every use check for damage of the applicators microdermabrasion head, defects can cause skin irritation or scratches and must be replaced.

How to use

Choose a treatment head (sensitive or normal). Place device tip on skin so as to feel the suction. Stretch the skin lightly with the fingers and move the device repeatedly 2-4 passes initially, lightly over the skin, continue for all areas to be treated. Ensure sensitive areas above and below the eyes are avoided. A treatment is approx. 5 minutes.

After the treatment

After the treatment your skin can accept a deeper penetration of skin care products. Do not forget to apply a sun protector the following 3 days before your skin is exposed to the sun.

Cleaning and care

Turn off device.
Clean devices treatment head after every use.
Device or charger must not be immersed in any fluid.
Do not use any type of abrasive pads or any aggressive cleaning fluids or paste.
Remove the treatment head from the device. This is magnetic and is easily done.
Clean under warm running water. With a small amount of dish washing liquid and a small brush the dirt particles are easily removed.
Can be left to dry it or with a soft microphase cloth.
Main device body and charger should also be wiped clean using damp microphase cloth.

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WM BEAUTYSYSTEMS AG & Co. KG